ENZYBREW L

Enzyme detergent for soaking in microbreweries

Enzyme detergent specially designed for cleaning bottling machines and removable parts (valves, rubber hoses, fittings,...) in microbreweries. ENZYBREW L breaks organic soiling such as protein waste and residues of polysaccharides completely down. With its phosphorus-free formula the product provides a deep and eco-friendly cleaning.

Advantages

- **Reinforced hygiene**: deep cleaning due to enzymes (removal of organic waste) and effective preparation for disinfection
- Removal of contamination risk
- Very easy to rinse: non-foaming product
- No danger for the user: non-corrosive (no risk of burns) and non-classified product
- **Guaranteed organoleptic quality** of finished product
- Environmental friendly alternative to caustic soda

Fields of application

Suitable for:

• Cleaning bottling machines and soaking removable parts (valves, rubber hoses, fittings, ...)







Directions for use

SOAKING

Suitable for cleaning after each use. Recommended dosage: 1%

- 1. Dilute the product in hot water (45-55°C) and mix well.
- 2. Soak removable parts for between 4 and 12 hours.
- 3. Rinse with hot water.
- 4. Disinfect.
- 5. Rinse with clean water.

Remark: If there is prolonged contact with light metals (aluminium, iron, copper), testing first on a small surface is recommended.

Temperature 45-55°C

Action soaking + rinsing

Concentration 1%

Contact time Between 4 and 12 hours

PH value

|--|

In solution

Certifications and regulations

Product compliant with:

- European Regulation EC 648/2004 concerning detergents
- Regulation EC 1272/2008 concerning the classification, packaging and labelling of substances and mixtures
- REACH European Regulation (EC 1907/2006) adopted to protect human health and the environment better from the risks linked to chemicals

Safety

For professional use only. Consult the safety sheet (available on <u>www.realco.be/msds</u>) and the product label for more information.

Storage

Storage in closed original package between +4°C and +25°C. Use the product within 12 months from the date of manufacture.

Packaging

Reference: LIQ1665 Box of 4x51 Pallet of 24 boxes (= 96 bottles)

